



Hot-start AptaTaq DNA polymerase

Cat: C101011 (250 U, 50 µl)

Cat: C101012 (500 U, 100 µl)

Store at -20 °C

Contents:

| Component | C101011 | C101012 |
|-----------------------------------|---------|---------|
| Apta Taq DNA poly. 5 U/µl | 250 U | 500 U |
| MgCl ₂ Solution 25 mM | 0.5 ml | 1 ml |
| 10X Buffer MgCl ₂ free | 0.5 ml | 1 ml |

Description:

This DNA polymerase is a mixture of Taq DNA polymerase and a temperature sensitive, aptamer-based inhibitor. The inhibitor binds reversibly to the enzyme, inhibiting polymerase activity at temperatures below 40 °C, but releases the enzyme during normal PCR cycling conditions. The aptamer-based hot start mechanism does not require a separate high temperature incubation step to activate the enzyme. The enzyme is inactive at room temperature, avoiding extension of non-specifically annealed primers or primer dimers and providing higher specificity of DNA amplification.

The activated enzyme maintains the same functionality as Taq DNA polymerase: it catalyzes 5' → 3' synthesis of DNA, has no detectable 3' → 5' proofreading exonuclease activity.

Kit storage:

This kit should be stored at -20 °C. Under this condition reagents are stable for two years from the date of production.

General Reaction Protocol:

1. Thaw 10X reaction buffer, dNTP mixture.
2. Mix the master mix thoroughly and dispense appropriate volumes into PCR tubes or plates.

| Component | Volume | Final Conc. |
|----------------------------------|----------|---------------|
| 10X Reaction Buffer | 2 µL | 1X |
| MgCl ₂ Solution 25 mM | 1.2 µL | 1.5 mM |
| 40 mM dNTPs Mix (10 mM each) | 0.5 µL | 0.25 mM |
| Upstream Primer (10 pmol/ µL) | 1 µL | 0.5 pmoles/µL |
| Downstream Primer (10 pmol/ µL) | 1 µL | 0.5 pmoles/µL |
| Template DNA | Variable | 10 fg~1 µg |
| PCR grade water | Variable | - |
| Apta Taq DNA poly. 5U/µl | 0.25 µL | 1.25 U |
| Total Volume | 20 µL | - |

3. Add templates DNA to the individual PCR tubes or wells containing the master mix.
4. Program the PCR machine according to the program outlined.

| Cycle | Time | Temp °C |
|-------|-----------|---------|
| 1 | 4 min | 95 |
| 30-35 | 30 sec | 94 |
| | 30 sec | 57 |
| | 30-60 sec | 72 |
| 1 | 5 min | 72 |

Notes:

Extension temperature is between 68 and 72°C. We highly recommend 68 °C for more efficiency of Pars Tous Taq DNA polymerase.

* For PCR products longer than 3~4 Kb, use an extension time of approximately 1 min per Kb DNA.

* A DNA fragment which is amplified by Taq DNA polymerase has A overhang, and it enables you to do cloning by using T-vector

Agarose gel Electrophoresis:

Run the total 5-7 µL of PCR products alongside 3µL DNA marker on a 2% agarose gel containing Green viewer DNA safe stain.

Disclaimers and Addresses:

This product is for **Research Use Only** and should only be used by trained professionals.

Unit 18, North East Food tech park, Tous Industrial Zone, Mashhad- IRAN
 Tel: (+98-51)35420843-5
 www.parstous.com

Fax: (+98-51) 35420846
 info@parstous.com