

Easy™ cDNA Synthesis Kit

Cat: A101161 (50 Reactions)
Cat: A101162 (100 Reactions)

Store at -20 °C

Contents:

Component	A101161	A101162
Buffer-Mix (2x)	500 μL	1000 μL
Enzyme Mix	100 μL	200 μL
DEPC-treated water	500 μL	1000 μL
cDNA Control Primer	50 μL	100 μL
Mix (B2M)		

Description:

Easy cDNA Synthesis kit contains all necessary components for conversion of total RNA or mRNA to **the single stranded cDNA**. The 2X Buffer mix solutions contains, RT buffer, 1mM dNTP mixture, 8mM MgCl₂, Oligo d(t)16, Random hexamer and stabilizer. Enzyme mix contains thermostable Hminus MMLV. RNase Inhibitor and stabilizer.

Features:

- Easy protocol & Minimum pipetting steps
- RNase minus MMLV enzyme

- Long mRNA synthesis
- High temperature reaction to destabilize RNA secondary structures

General Reaction Protocol (first strand cDNA synthesis):

1- Mix the template RNA (total RNA or Poly (A) mRNA) and other kit components in RNase-free tube as below table.

DEPC-treated water	Up to 20 µL
Enzyme Mix	2
Buffer-Mix (2x)	10
Template RNA 1 ng~5 μg	ХμL

- 2- Mix the above mixture by quick vortex.
- 3- Incubate 10 min at 25 °C.
- 4. Incubate 60 min at 47 °C.
- 5. Stop the reaction by heating at 85 $^{\circ}$ C for 5 minutes. Chill on the ice or at 4 $^{\circ}$ C.

Note: To perform PCR, you can add the finished RT reaction up to 1:5 of the final PCR volume.

Component	Vol. (μL)
Taq 2X Premix	10
cDNA Control primer mix	1
cDNA	1-5
PCR Grade Water	Up to 20 μL

cDNA Control PCR Reaction:

- 1- Prepare a reaction mix according to the table.
- 2- For negative tube use PCR grade water. The final volume in each PCR reaction tubes is 20 µl.

Cycle	Time	Temp °C
1	4 min	95
·	30 sec	94
35	30 sec	57
	30 sec	72
1	5 min	72

Note: It is recommended that all of the PCR components be premixed in a sufficient quantity for daily needs and then dispensed into the individual reaction tubes.

Amplification protocol:

Agarose gel Electrophoresis:

Run the total 5-7 μ L of PCR products alongside 3 μ L DNA marker on a 2% agarose gel containing Green viewer DNA safe stain.

The B2M primers amplify a band of approximately 230 bp from human. mouse and rat B2M cDNA.

Disclaimers and Addresses:

This product is for **Research Use Only** and should only be used by trained professionals.

Unit 18, North East Food tech park, Tous Industrial Zone, Mashhad-IRAN

Tel: (+98-51)35420843-5 www.parstous.com Fax: (+98-51) 35420846 info@parstous.com